

Experiment Results Report

Experiment Summary

Experiment Name: CoV2 400-058 HMP QuantS_v4
Experiment Type: Standard Curve
Chemistry: TaqMan® Reagents
BarCode:
File Name: CoV2 400-058 HMP QuantS_v4.edt
Run Started: Not Started
Run Finished:
Run Duration:
Date EDT Modified: 06-24-2020 18:22:48 CDT
Date EDT Created: 04-22-2020 09:59:14 CDT
User Name:
Number of Wells Used: 0
Number of Wells with Results: 0
Instrument Name:
Instrument Type: QuantStudio™ 5 System
Instrument Serial Number:
Model/Block Type: QuantStudio™ 5 System / 96-Well 0.2-mL Block
Quantification Cycle Setting: Ct
Stage/Step for Analysis: Stage2, Step3
Comments:



Experiment Results Report

Reagent Information

Type	Name	Part Number	Lot Number	Expiration Date
------	------	-------------	------------	-----------------



For Research Use Only, not for use in diagnostic procedures.

Experiment Results Report

TC Protocol

Block Type: 96-Well 0.2-mL Block
 Run Mode: Standard
 Sample Volume: 15.0
 Cover Temperature: 105.0
 PCR Filter Set: m1,x1 m2,x1 m2,x2 m3,x2 m3,x3 m4,x3 m4,x4 m5,x4 m5,x5 m6,x5 m6,x6
 Melt Curve Filter Set: Not Applicable

Stage	No. of Repetitions	Starting Cycle	Auto Delta Enabled
Hold Stage	1	1	false
Cycling Stage	40	1	false

Stage	Step	Collection Flag	Ramp Rate	Temperature	Hold Time	Auto Delta Temperature	Auto Delta Hold Time
Hold Stage	1	Disabled	1.6 °C/s	43.0	900	0.0	0
Hold Stage	2	Disabled	1.6 °C/s	95.0	120	0.0	0
Cycling Stage	1	Disabled	1.6 °C/s	95.0	3	0.0	0
Cycling Stage	2	Disabled	1.6 °C/s	68.0	10	0.0	0
Cycling Stage	3	Enabled	1.6 °C/s	74.0	5	0.0	0

Customer is responsible for any validation of assays and compliance with regulatory requirements that pertain to their procedures and uses of reagents and/or instruments. ©2018 Thermo Fisher Scientific. All rights reserved. The trademarks mentioned herein are the property of Thermo Fisher Scientific or their respective owners.

