

Date/Time: 2/9/2017 1:51:03 AM  
 Username: ADMIN

# UDP Report

Instrument #: MX0000(1)  
 Version #: V4.50 A

## PCP-350010AMAX

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time

PCP-350010AMAX	4	ExK DNA-1 (Urine)	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing	20
		Alias	PCP	--	--	SPC	62	1600
		Gain	60	--	--	60	700	500
		Threshold	100	--	--	100	--	--

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	120	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 45

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	46.7	62	Yes

Performed By: \_\_\_\_\_ Date: \_\_\_\_\_ Printed: 2017-02-09 01:51:03  
 Reviewed By: \_\_\_\_\_ Date: \_\_\_\_\_ Page: 1/1

