

Date/Time: 10/18/2018 5:14:32 PM  
 Username: ADMIN

# UDP Report

Instrument #: CT0156(1)  
 Version #: V4.70 A

## 400005CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time
							Liquid Level Sensing Enabled	
400005CMAX S4	1 ExK DNA-2 - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer					Call Ct at Threshold Crossing	9
--	Alias	Strep	Neis	Haem	--	SPC	62	1600
	Gain	60	60	60	--	70	937.5	500
	Threshold	250	250	250	--	250	--	--
							Enabled	

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	13.7	99	No
--	45	62	Yes

		Strep			
Strep - PCR	✓	✓	✗	✗	
SPC - PCR	✓	✗	✓	✗	
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR	

		Neis			
Neis - PCR	✓	✓	✗	✗	
SPC - PCR	✓	✗	✓	✗	
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR	

		Haem			
Haem - PCR	✓	✓	✗	✗	
SPC - PCR	✓	✗	✓	✗	
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR	



Performed By: \_\_\_\_\_

Date: \_\_\_\_\_

Printed: 2018-10-18 17:14:33

Reviewed By: \_\_\_\_\_

Date: \_\_\_\_\_

Page: 1/1