

Date/Time: **11/27/2018 7:48:27 AM**  
 Username: **ADMIN**

# UDP Report

Instrument #: **CT0156(1)**  
 Version #: **V4.70 A**

## 400014AMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation		Lysis Heat Time		
			Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
										Sample Volume	Wash Volume
										Neutralization Volume	DNase Heat Time
										Liquid Level Sensing Enabled	
400014AMAX S4	1	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing		20		
--		Alias	PCP	--	--	--	SPC	62	1600		
		Gain	60	--	--	--	60	700	500		
		Threshold	100	--	--	--	100	--	--		
								Enabled			

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 45

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	46.7	62	Yes

PCP				
PCP - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reviewed By: \_\_\_\_\_ Date: \_\_\_\_\_

Printed: 2018-11-27 07:48:28