

Date/Time: **10/23/2018 9:53:19 AM**  
 Username: **ADMIN**

# UDP Report

Instrument #: **CT0156(1)**  
 Version #: **V4.70 A**

## 400020CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation		Lysis Heat Time		
			Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
										Sample Volume	Wash Volume
										Neutralization Volume	DNase Heat Time
										Liquid Level Sensing Enabled	
400020CMAX S4	1	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing		30		
--		Alias	--	pIS1001	IS481	--	SPC	62	1600		
--		Gain	--	60	60	--	60	500	500		
--		Threshold	--	200	200	--	200	--	--		
								Enabled			

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	44.6	60	Yes

pIS1001				
pIS1001 - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

IS481				
IS481 - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: \_\_\_\_\_

Date: \_\_\_\_\_

Printed: 2018-10-23 09:53:20

Reviewed By: \_\_\_\_\_

Date: \_\_\_\_\_

Page: 1/1