

Date/Time: 10/25/2018 9:33:08 AM  
 Username: ADMIN

# UDP Report

Instrument #: CR0009(1)  
 Version #: V4.70 A

## 400028CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time	
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
							400028CMAX S4	1	ExK DNA-3 - 4-Snap
		Alias	Efaec	GAS	KlebOP	Ecoli	SPC	62	1600
		Gain	60	60	60	60	60	937.5	500
		Threshold	200	200	200	300	200	--	--
		Excitation Channel 475/520	--	0	0	0	0		
--		Excitation Channel 530/565	3	--	0	0	0		
		Excitation Channel 585/630	0	0	--	5	0		
		Excitation Channel 630/665	0	0	5	--	0		
		Excitation Channel 680/715	0	0	0	0	--		

Enabled

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	97	No

**Step Name:** New Step  
**Profile Type:** 3 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	12	97	No
--	24.6	60	Yes
--	10.7	72	No

Efaec				
Efaec - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
<b>Results</b>	➕ POS	➕ POS	➖ NEG	⚠ UNR

GAS				
GAS - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
<b>Results</b>	➕ POS	➕ POS	➖ NEG	⚠ UNR

KlebOP				
KlebOP - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
<b>Results</b>	➕ POS	➕ POS	➖ NEG	⚠ UNR

Ecoli				
Ecoli - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
<b>Results</b>	➕ POS	➕ POS	➖ NEG	⚠ UNR



Performed By: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reviewed By: \_\_\_\_\_ Date: \_\_\_\_\_

Printed: 2018-10-25 09:33:08