

Date/Time: **10/24/2018**
10:31:23 AM
 Username: **ADMIN**

UDP Report

Instrument #: **CT0156(1)**
 Version #: **V4.70 A**

400035CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation		Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
							Call Ct at Threshold Crossing		
400035CMAX S4	1	ExK DNA-3 - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer					9	
		Alias	Glab	Alb	Krus	TV	SPC	1600	
		Gain	60	80	60	60	60	500	
		Threshold	200	200	300	200	200	--	
		Excitation Channel 475/520	--	0	0	0	0	--	
--		Excitation Channel 530/565	5	--	5	0	0		
		Excitation Channel 585/630	0	0	--	5	0		
		Excitation Channel 630/665	0	0	5	--	0		
		Excitation Channel 680/715	0	0	0	3	--		

Enabled

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	14.2	99	No
--	47.2	65	Yes

	Glab			
Glab - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➔ NEG	⚠ UNR

	Alb			
Alb - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➔ NEG	⚠ UNR

	Krus			
Krus - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➔ NEG	⚠ UNR

	TV			
TV - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➔ NEG	⚠ UNR



Performed By: _____

Date: _____

Printed: 2018-10-24 10:31:23

Reviewed By: _____

Date: _____

Page: 1/1