

Date/Time: **10/18/2018**  
**11:28:29 AM**  
 Username: **ADMIN**

# UDP Report

Instrument #: **CT0156(1)**  
 Version #: **V4.70 A**

## 400007CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time
							Liquid Level Sensing Enabled	
400007CMAX S4	1 ExK DNA-2 - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer					Call Ct at Threshold Crossing	9
--	Alias	vanA	vanC1/2	vanB	--	SPC	62	1600
	Gain	60	40	80	--		937.5	500
	Threshold	100	400	100	--	100	--	--
							Enabled	

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	12.8	99	No
--	46	62	Yes

vanA				
vanA - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

vanC1/2				
vanC1/2 - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

vanB				
vanB - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: \_\_\_\_\_

Date: \_\_\_\_\_

Printed: 2018-10-18 11:28:30

Reviewed By: \_\_\_\_\_

Date: \_\_\_\_\_

Page: 1/1