

Date/Time: **10/22/2018**
12:32:20 PM
 Username: **ADMIN**

UDP Report

Instrument #: **CT0156(1)**
 Version #: **V4.70 A**

400030CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time	
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
400030CMAX S4	1	ExK DNA-1 (Plasma/Serum) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing	10	
--		Alias	--	--	EBV	--	SPC	37	1600
--		Gain	--	--	60	--	60	700	500
--		Threshold	--	--	200	--	200	--	--
							Enabled		

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	14.3	99	No
--	44.5	62	Yes

EBV				
EBV - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: _____

Date: _____

Printed: 2018-10-22 12:32:20

Reviewed By: _____

Date: _____

Page: 1/1