

Date/Time: 10/25/2018 9:37:43 AM
 Username: ADMIN

UDP Report

Instrument #: CT0271(2)
 Version #: V4.70 A

400038CMAX S4

| Test Name | Version | Extraction Type | Master Mix Format | | | | Ct Calculation | Lysis Heat Time | |
|---------------|------------------|--------------------|--|---------|---------|---------|-------------------------------|------------------------------|-----------------|
| Notes | Channel Settings | 475/520 | 530/565 | 585/630 | 630/665 | 680/715 | Lysis Temperature | Sample Tip Height | |
| | | | | | | | | Sample Volume | Wash Volume |
| | | | | | | | | Neutralization Volume | DNase Heat Time |
| | | | | | | | | Liquid Level Sensing Enabled | |
| 400038CMAX S4 | 1 | ExK DNA-3 - 4-Snap | Type 4: Lyophilized MM with Rehydration Buffer | | | | Call Ct at Threshold Crossing | 9 | |
| -- | | Alias | -- | -- | Tpal | -- | SPC | 62 | 1600 |
| | | Gain | -- | -- | 60 | -- | 60 | 500 | 500 |
| | | Threshold | -- | -- | 200 | -- | 200 | -- | -- |
| | | | | | | | | Enabled | |

Step Name: New Step
Profile Type: Hold **Cycles:** 1

| Type | Time (s) | Temp (°C) | Detect |
|------|----------|-----------|--------|
| -- | 300 | 99 | No |

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

| Type | Time (s) | Temp (°C) | Detect |
|------|----------|-----------|--------|
| -- | 13.7 | 99 | No |
| -- | 45 | 62 | Yes |

| Tpal | | | | |
|------------|-------|-------|-------|-------|
| Tpal - PCR | ✓ | ✓ | ✗ | ✗ |
| SPC - PCR | ✓ | ✗ | ✓ | ✗ |
| Results | ⊕ POS | ⊕ POS | ⊖ NEG | ⚠ UNR |



Performed By: _____

Date: _____

Printed: 2018-10-25 09:37:43

Reviewed By: _____

Date: _____

Page: 1/1