

BioGX
Molecular Made Easy



REAGENT CATALOG

Syndromic Multiplexes (RUO)



Syndromic Group	Product Name	Targets	Reactions per kit	Part No.
Antimicrobial Resistance	vanA, vanB, vanC1, vanC2/3	vanA, vanB, vanC1, vanC2/3, SPC	24	450-016-C-MAX
	vanA, vanB	vanA, vanB, SPC	24	450-017-C-MAX
	vanA, vanB	vanA, vanB, RNase P	64	450-017-LMP
	KPC, NDM, OXA-48	KPC, NDM, OXA-48, SPC	24	450-018-C-MAX
	ABRM1 - KPC, NDM, OXA-48	KPC, NDM, OXA-48, RNase P	64	450-018-LMP
	VIM, GES, mcr-1	VIM, GES, mcr-1, SPC	24	450-019-C-MAX
	ABRM2 - VIM, GES, mcr-1	VIM, GES, mcr-1, RNase P	64	450-019-LMP
	ABRM3 - vanA, vanB, mecA	vanA, vanB, mecA, RNase P	64	450-081-LMP
	ABRS1 - Aminoglycoside Resistance (Aad and AAC)	Aad and AAC	64	450-098-PNP
	ABRS2 - Aminoglycoside Resistance (APH)	APH	64	450-099-PNP
	ABRS3 - Cephalosporin Resistance (ACC, ACT, FOX, MIR)	ACC, ACT, FOX, MIR	64	450-100-PNP
	ABRS4 - Cephalosporin Resistance (CMY/DHA/LAT/MOX/MOR and BIL)	CMY/DHA/LAT/MOX/MOR and BIL	64	450-101-PNP
	ABRS5 - ESBL (CTX-M)	CTX-M	64	450-102-PNP
	ABRS6 - ESBL (SHV)	SHV	64	450-103-PNP
	ABRS7 - ESBL (TEM)	TEM	64	450-104-PNP
	ABRS8 - Carbapenem Resistance (IMP)	IMP	64	450-105-PNP
	ABRS9 - Trimethoprim Resistance (dfrA1/dfrA5/dfrA12/dfrA17)	dfrA1/dfrA5/dfrA12/dfrA17	64	450-106-PNP
ABRS10 - Sulfonamide Resistance (sul1/sul2)	sul1/sul2	64	450-107-PNP	
ABRS11 - Nitroimidazole Resistance (nimBDEJ)	nimBDEJ	64	450-108-PNP	
ABRS12 - Macrolide Resistance (ermABC)	ermABC	64	450-109-PNP	
ABRS13 - Quinolone Resistance (qnrAS)	qnrAS	64	450-110-PNP	
ABRS14 - Quinolone Resistance (qnrB)	qnrB	64	450-111-PNP	
ABRS15 - Tetracycline Resistance (tetABMS)	tetABMS	64	450-112-PNP	
SPC - RNase P	RNase P	64	450-113-PNP	
Bloodborne	Plasmodium spp., P. falciparum, P. vivax	Plasmodium spp., P. falciparum, P. vivax, RNase P	24	450-032-D-MAX
	Plasmodium spp., P. falciparum, P. vivax	Plasmodium spp., P. falciparum, P. vivax, RNase P	64	450-032-LMP
	Anaplasma phagocytophilum, Babesia microti, Ehrlichia spp.	Anaplasma phagocytophilum, Babesia microti, Ehrlichia spp., SPC	24	450-033-D-MAX
COVID-19	COVID-19, Flu A, Flu B, RSV	COVID-19, Flu A, Flu B, RSV, SPC	24	450-015-G-MAX
	COVID-19, Flu A, Flu B, RSV	COVID-19, Flu A, Flu B, RSV, IAC	64	450-015-H-LMP
	Xfree COVID-19, Flu A, Flu B	SARS CoV-2 (N1 and RdRp genes), Flu A, Flu B, IAC	16	450-090-H-PXL
Gastrointestinal	Campylobacter spp., Shigella spp., Salmonella spp.	Campylobacter spp., Shigella spp., Salmonella spp., SPC	24	450-020-C-MAX
	Escherichia coli O157, stx-1, stx-2	Escherichia coli O157, stx-1, stx-2, SPC	24	450-021-C-MAX
	Cryptosporidium parvum, Giardia intestinalis, Entamoeba histolytica, Dientamoeba fragilis	Cryptosporidium parvum, Giardia intestinalis, Entamoeba histolytica, Dientamoeba fragilis, SPC	24	450-022-C-MAX
	Cyclospora cayetenensis, Microsporidia (Encephalitozoon cuniculi, Encephalitozoon hellem, Encephalitozoon intestinalis, Enterocytozoon bieneusi)	Cyclospora cayetenensis, Microsporidia (Encephalitozoon cuniculi, Encephalitozoon hellem, Encephalitozoon intestinalis, Enterocytozoon bieneusi), SPC	24	450-023-C-MAX
	Microsporidia (Encephalitozoon cuniculi, Encephalitozoon hellem, Encephalitozoon intestinalis, Enterocytozoon bieneusi)	Microsporidia (Encephalitozoon cuniculi, Encephalitozoon hellem, Encephalitozoon intestinalis, Enterocytozoon bieneusi), RNase P	64	450-063-LMP
Immunocompromised, Meningitis, Encephalitis	Enterovirus	Enterovirus, SPC	24	450-006-C-MAX
	Enterovirus	Enterovirus, RNase P, IAC	64	450-006-LMP
	Neisseria meningitidis, Streptococcus pneumoniae, Haemophilus influenzae	Neisseria meningitidis, Streptococcus pneumoniae, Haemophilus influenzae, SPC	24	450-024-C-MAX
	Neisseria meningitidis, Streptococcus pneumoniae, Haemophilus influenzae	Neisseria meningitidis, Streptococcus pneumoniae, Haemophilus influenzae, RNase P	64	450-024-LMP
	Escherichia coli, Listeria spp., Streptococcus agalactiae	Escherichia coli, Listeria spp., Streptococcus agalactiae, SPC	24	450-025-C-MAX
	Escherichia coli, Listeria spp., Streptococcus agalactiae	Escherichia coli, Listeria spp., Streptococcus agalactiae, RNase P	64	450-025-LMP
	BK Virus	BK Virus, SPC	24	450-026-C-MAX
	CMV	CMV, SPC	24	450-027-C-MAX
	JC Virus	JC Virus, SPC	24	450-028-C-MAX
	Epstein-Barr Virus	Epstein-Barr Virus, SPC	24	450-029-C-MAX
HSV-1, HSV-2, VZV	HSV-1, HSV-2, VZV, SPC	24	450-030-C-MAX	

Research Use Only Products (450-series part numbers). Not for Use in Diagnostic Procedures.

Controls and other info available at <https://www.biogx.com/syndromic-multiplexes/>

Syndromic Multiplexes (RUO)

Syndromic Group	Product Name	Targets	Reactions per kit	Part No.
Immunocompromised, Meningitis, Encephalitis	Xfree HSV-1, HSV-2, VZV	HSV-1, HSV-2, VZV, RNase P	16	450-030-E-PXL
	HSV-1, HSV-2, VZV	HSV-1, HSV-2, VZV, RNase P	64	450-030-LMP
	HSV-1, HSV-2, VZV, CMV	HSV-1, HSV-2, VZV, CMV, SPC	24	450-031-C-MAX
	Xfree hMPXV	Monkeypox Virus, Orthopox Virus, RNase P	16	450-086-E-PXL
	Acanthamoeba spp., Balamuthia mandrillaris, Naegleria fowleri	Acanthamoeba spp., Balamuthia mandrillaris, Naegleria fowleri, SPC	24	450-095-C-MAX
	Acanthamoeba spp., Balamuthia mandrillaris, Naegleria fowleri	Acanthamoeba spp., Balamuthia mandrillaris, Naegleria fowleri, SPC	24	450-095-C-MAX
Nail Fungus	Candida albicans, Candida glabrata, Candida tropicalis/ Candida parapsilosis, Candida krusei	Candida albicans, Candida glabrata, Candida tropicalis/ Candida parapsilosis, Candida krusei, SPC	24	450-036-C-MAX
	Dermatophyte spp., Fungus spp., Trichophyton rubrum	Dermatophyte spp., Fungus spp., Trichophyton rubrum, SPC	24	450-048-C-MAX
Respiratory	Adenovirus (species A-G), Bordetella pertussis, Bordetella parapertussis	Adenovirus (species A-G), Bordetella pertussis, Bordetella parapertussis, SPC	24	450-002-C-MAX
	Bordetella pertussis, B. parapertussis, B. holmesii, pertussis toxin	Bordetella pertussis, B. parapertussis, B. holmesii, pertussis toxin, SPC	24	450-003-C-MAX
	Bordetella pertussis, B. parapertussis	Bordetella pertussis, B. parapertussis, SPC	24	450-004-C-MAX
	Coccidioides spp.	Coccidioides spp., SPC	24	450-005-W-MAX
	Enterovirus	Enterovirus, SPC	24	450-006-C-MAX
	Enterovirus	Enterovirus, RNase P, IAC	64	450-006-LMP
	Mycoplasma pneumoniae, Legionella pneumophila, Chlamydia pneumoniae	Mycoplasma pneumoniae, Legionella pneumophila, Chlamydia pneumoniae, SPC	24	450-007-A-MAX
	Mycoplasma pneumoniae, Legionella spp., Chlamydia pneumoniae, Chlamydia psittaci	Mycoplasma pneumoniae, Legionella spp., Chlamydia pneumoniae, Chlamydia psittaci, SPC	24	450-008-C-MAX
	Mycoplasma pneumoniae, Legionella pneumophila, Chlamydia pneumoniae, Pneumocystis jirovecii	Mycoplasma pneumoniae, Legionella pneumophila, Chlamydia pneumoniae, Pneumocystis jirovecii, SPC	24	450-009-C-MAX
	Mycobacterium tuberculosis Complex	Mycobacterium tuberculosis Complex, SPC [MTB-M. africanum, M. bovis, M. canettii, M. caprae, M. microti, M. orygis, M. pinnipedii and M.tuberculosis]	24	450-010-W-MAX
	Aspergillus spp.	Aspergillus spp., SPC	24	450-011-C-MAX
	Pneumocystis jirovecii	Pneumocystis jirovecii, SPC	24	450-012-A-MAX
	Streptococcus pyogenes	Streptococcus pyogenes, SPC	24	450-013-C-MAX
	Streptococcus pyogenes	Streptococcus pyogenes, RNase P	64	450-013-LMP
	Flu A, Flu B, RSV A/B	Flu A, Flu B, RSV A/B, SPC	24	450-014-C-MAX
	Mycobacterium tuberculosis complex, Mycobacterium avium complex, Mycobacterium spp.	Mycobacterium tuberculosis complex, Mycobacterium avium complex, Mycobacterium spp., SPC	24	450-054-C-MAX
	Mycobacterium tuberculosis complex, Mycobacterium avium complex, Mycobacterium spp.	Mycobacterium tuberculosis complex, Mycobacterium avium complex, Mycobacterium spp., RNase P	64	450-054-LMP
	Mycobacterium chelonae, Mycobacterium abscessus group, Mycobacterium fortuitum complex	Mycobacterium chelonae, Mycobacterium abscessus group, Mycobacterium fortuitum complex, SPC	24	450-055-C-MAX
	Mycobacterium chelonae, Mycobacterium abscessus group, Mycobacterium fortuitum complex	Mycobacterium chelonae, Mycobacterium abscessus group, Mycobacterium fortuitum complex, RNase P	64	450-055-LMP
	Flu A, Flu B, RSV A/B	Flu A, Flu B, RSV A/B, IAC	150	450-070-HMP
Xfree Group A Streptococcus	Streptococcus pyogenes, RNase P	16	450-092-PXL	
Xfree Mycoplasma pneumoniae	Mycoplasma pneumoniae, RNase P	16	450-093-E-PXL	
Legionella pneumophila, Legionella spp.	Legionella pneumophila, Legionella spp., SPC	24	450-094-C-MAX	
Legionella pneumophila, Legionella spp.	Legionella pneumophila, Legionella spp., RNase P	64	450-094-LMP	
Sexual Health	Xfree HSV-1, HSV-2, VZV	HSV-1, HSV-2, VZV, RNase P	16	450-030-E-PXL
	Candida albicans, Candida glabrata, Candida krusei, Trichomonas vaginalis	Candida albicans, Candida glabrata, Candida krusei, Trichomonas vaginalis, SPC	24	450-035-C-MAX
	Candida albicans, Candida glabrata, Candida krusei, Trichomonas vaginalis	Candida albicans, Candida glabrata, Candida krusei, Trichomonas vaginalis, RNase P	64	450-035-LMP
	Gardnerella vaginalis, Trichomonas vaginalis, Candida albicans, Candida glabrata	Gardnerella vaginalis, Trichomonas vaginalis, Candida albicans, Candida glabrata, SPC	24	450-037-C-MAX
	HSV-1, HSV-2, Haemophilus ducreyi	HSV-1, HSV-2, Haemophilus ducreyi, SPC	24	450-038-C-MAX
	Mycoplasma genitalium, M. hominis, Ureaplasma urealyticum, U. parvum	Mycoplasma genitalium, M. hominis, Ureaplasma urealyticum, U. parvum, SPC	24	450-039-C-MAX
	Treponema pallidum	Treponema pallidum, SPC	24	450-040-C-MAX
	Xfree Treponema pallidum	Treponema pallidum, RNase P	16	450-040-E-PXL
	Treponema pallidum	Treponema pallidum, RNase P	64	450-040-LMP

Research Use Only Products (450-series part numbers). Not for Use in Diagnostic Procedures.

Controls and other info available at <https://www.biogx.com/syndromic-multiplexes/>

Syndromic Multiplexes (RUO)



Syndromic Group	Product Name	Targets	Reactions per kit	Part No.
Sexual Health	Gardnerella vaginalis, Trichomonas vaginalis, Candida spp.	Gardnerella vaginalis, Trichomonas vaginalis, Candida spp., RNase P	64	450-078-LMP
	SRM10 - HSV-1, HSV-2	HSV-1, HSV-2, RNase P	64	450-079-LMP
	Chlamydia trachomatis, Neisseria gonorrhoeae, Enterococcus faecalis	Chlamydia trachomatis, Neisseria gonorrhoeae, Enterococcus faecalis, RNase P	64	450-080-LMP
	Xfree hMPXV	hMPXV, RNase P	16	450-086-E-PXL
	SRM13 - Mycoplasma genitalium, Mycoplasma hominis, Ureaplasma spp.	Mycoplasma genitalium, Mycoplasma hominis, Ureaplasma spp., RNase P	64	450-087-LMP
	SRM14 - Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis	Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis, RNase P	64	450-088-LMP
Urinary Tract/Wound	Acinetobacter baumannii, Enterobacter spp., Morganella morganii, Kingella kingae	Acinetobacter baumannii, Enterobacter spp., Morganella morganii, Kingella kingae, SPC	24	450-041-C-MAX
	Bacteroides fragilis, Citrobacter freundii, Clostridium novyi/Clostridium septicum, Clostridium perfringens	Bacteroides fragilis, Citrobacter freundii, Clostridium novyi/Clostridium septicum, Clostridium perfringens, SPC	24	450-042-C-MAX
	Candida auris	Candida auris, SPC	24	450-043-C-MAX
	Candida auris	Candida auris, RNase P	64	450-043-LMP
	Enterococcus faecalis, Enterococcus faecium, Streptococcus pyogenes, Streptococcus agalactiae	Enterococcus faecalis, Enterococcus faecium, Streptococcus pyogenes, Streptococcus agalactiae, SPC	24	450-044-C-MAX
	Escherichia coli, Klebsiella spp., Staphylococcus saprophyticus, Enterococcus spp.	Escherichia coli, Klebsiella spp., Staphylococcus saprophyticus, Enterococcus spp., SPC	24	450-045-C-MAX
	Escherichia coli, Klebsiella spp., Pseudomonas aeruginosa, Proteus spp.	Escherichia coli, Klebsiella spp., Pseudomonas aeruginosa, Proteus spp., SPC	24	450-046-C-MAX
	Escherichia coli, Klebsiella spp, Enterococcus faecalis, Streptococcus pyogenes	Escherichia coli, Klebsiella spp, Enterococcus faecalis, Streptococcus pyogenes, SPC	24	450-047-C-MAX
	Staphylococcus aureus, mecA, mecC, Streptococcus agalactiae	Staphylococcus aureus, mecA, mecC, Streptococcus agalactiae, SPC	24	450-049-C-MAX
	Staphylococcus epidermidis group, Staphylococcus aureus, Staphylococcus lugdunensis, mecA	Staphylococcus epidermidis group, Staphylococcus aureus, Staphylococcus lugdunensis, mecA, SPC [S. epidermidis group - S. hominis, S. simulans, S. warneri, S. epidermidis, S. haemolyticus]	24	450-050-C-MAX
	SRM1 - Bacteroides fragilis, Clostridium novyi/Clostridium septicum, Clostridium perfringens	Bacteroides fragilis, Clostridium novyi/Clostridium septicum, Clostridium perfringens, RNase P	64	450-068-LMP
	SRM2 - Staphylococcus aureus, mecA, mecC	Staphylococcus aureus, mecA, mecC, RNase P	64	450-069-LMP
	SRM3 - Acinetobacter baumannii, Morganella morganii, Escherichia coli	Acinetobacter baumannii, Morganella morganii, Escherichia coli, RNase P	64	450-071-LMP
	SRM4 - Klebsiella aerogenes, Enterobacter cloacae complex, Proteus vulgaris	Klebsiella aerogenes, Enterobacter cloacae complex, Proteus vulgaris, RNase P	64	450-072-LMP
	SRM5 - Citrobacter freundii, Citrobacter koseri, Salmonella sp.	Citrobacter freundii, Citrobacter koseri, Salmonella sp., RNase P	64	450-073-LMP
	SRM6 - Klebsiella oxytoca, Klebsiella pneumoniae, Pseudomonas aeruginosa	Klebsiella oxytoca, Klebsiella pneumoniae, Pseudomonas aeruginosa, RNase P	64	450-074-LMP
SRM7 - Proteus mirabilis, Providencia stuartii, Serratia marcescens	Proteus mirabilis, Providencia stuartii, Serratia marcescens, RNase P	64	450-075-LMP	
SRM8 - Staphylococcus aureus, Staphylococcus lugdunensis, Staphylococcus saprophyticus	Staphylococcus aureus, Staphylococcus lugdunensis, Staphylococcus saprophyticus, RNase P	64	450-076-LMP	
SRM9 - Staphylococcus epidermidis, Staphylococcus haemolyticus, Streptococcus agalactiae	Staphylococcus epidermidis, Staphylococcus haemolyticus, Streptococcus agalactiae, RNase P	64	450-077-LMP	
SRM10 - HSV-1, HSV-2	HSV-1, HSV-2, RNase P	64	450-079-LMP	
SRM11 - Enterococcus faecium, Enterococcus faecalis	Enterococcus faecium, Enterococcus faecalis, RNase P	64	450-082-LMP	
SRM12 - Gardnerella vaginalis, Candida spp.	Gardnerella vaginalis, Candida spp., RNase P	64	450-083-LMP	
SRM13 - Mycoplasma genitalium, Mycoplasma hominis, Ureaplasma spp.	Mycoplasma genitalium, Mycoplasma hominis, Ureaplasma spp., RNase P	64	450-087-LMP	
SRM14 - Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis	Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis, RNase P	64	450-088-LMP	
SRM15 - Candida glabrata, Candida auris, Candida tropicalis	Candida glabrata, Candida auris, Candida tropicalis, RNase P	64	450-096-LMP	
SRM16 - Candida krusei, Candida albicans, Candida parapsilosis	Candida krusei, Candida albicans, Candida parapsilosis, RNase P	64	450-097-LMP	
Women's Health	Atopobium vaginae, Mobiluncus curtisii, Mobiluncus mulieris, Prevotella bivia	Atopobium vaginae, Mobiluncus curtisii, Mobiluncus mulieris, Prevotella bivia, SPC	24	450-034-C-MAX
	Candida albicans, Candida glabrata, Candida tropicalis/ Candida parapsilosis, Candida krusei	Candida albicans, Candida glabrata, Candida tropicalis/ Candida parapsilosis, Candida krusei, SPC	24	450-036-C-MAX
Sample-Ready™ Master Mixes	Sample-Ready DNA Master Mix Lyophilized Open System Reagents	not applicable	16 (direct) / 64 (extracted)	450-051-XMM
	Sample-Ready RNA Master Mix Lyophilized Open System Reagents	not applicable	16 (direct) / 64 (extracted)	450-058-XMM

Syndromic Multiplexes (CE-IVD)

Syndromic Group	Product Name	Targets	No. of Reactions per kit	Part No.
Antimicrobial Resistance	Vancomycin Resistance	vanA, vanB, vanC1, vanC2/3, SPC	24	400-007-C-MAX
	Carbapenem Resistance KNO	KPC, NDM-1, OXA-48, SPC	24	400-015-C-MAX
	Carbapenem/Colistin Resistance VGM	VIM, GES, mcr-1, SPC	24	400-017-C-MAX
	vanA/B	vanA, vanB, SPC	24	400-037-C-MAX
Gastrointestinal	Parasitic Enteric CGED	Cryptosporidium parvum, Giardia intestinalis, Entamoeba histolytica, Dientamoeba fragilis, SPC	24	400-009-C-MAX
Immunocompromised, Meningitis, Encephalitis	Bacterial Meningitis NSH	Neisseria meningitidis, Streptococcus pneumoniae, and Haemophilus influenzae, SPC	24	400-005-C-MAX
	Bacterial Meningitis ELGBS	Escherichia coli, Listeria spp., Group B Streptococcus, SPC	24	400-006-C-MAX
	EBV	Epstein-Barr Virus, SPC	24	400-030-C-MAX
	Viral Meningitis HSV-1, HSV-2, VZV	Herpes Simplex Virus 1, Herpes Simplex Virus 2, Varicella Zoster Virus, SPC	24	400-057-C-MAX
Respiratory	Bordetella Specification plus toxin	Bordetella pertussis, Bordetella parapertussis, Bordetella holmesii, pertussis toxin, SPC	24	400-002-C-MAX
	Mycobacterium tuberculosis Complex	Mycobacterium tuberculosis Complex, SPC	24	400-010-W-MAX
	Pneumocystis jirovecii	Pneumocystis jirovecii pneumonia, SPC	24	400-014-A-MAX
	Bordetella pertussis/parapertussis	Bordetella pertussis, Bordetella parapertussis, SPC	24	400-020-C-MAX
	Group A Streptococcus	Streptococcus pyogenes, SPC	24	400-033-C-MAX
	Flu A, Flu B, RSV A/B	Influenza A, Influenza B, Respiratory Syncytial Virus A/B, RNase P	24	400-056-C-MAX
	SARS-CoV-2 HMP - High Volume Multi-Platform PCR Reagents	SARS CoV-2 N1, SARS CoV-2 N2, RNase P)	24 (BD MAX) 48 (CFX96, QS5)	400-058-E-HMP
	COVID-19, Flu A, Flu B, RSV	SARS CoV-2 N1, Influenza A, Influenza B, Respiratory Syncytial Virus A/B	24	400-060-G-MAX
Sexual Health	Mycoplasma Ureaplasma	Mycoplasma genitalium, Ureaplasma urealyticum, Mycoplasma hominis, Ureaplasma parvum, SPC	24	400-003-C-MAX
	Candida Specification GATPK	Candida glabrata, Candida albicans, Candida tropicalis, Candida parapsilosis, Candida krusei, SPC	24	400-012-C-MAX
	Lesion HSVHD	Herpes Simplex Virus 1, Herpes Simplex Virus 2, Haemophilus ducreyi, SPC	24	400-027-C-MAX
	Candida-Trichomonas	Candida albicans, Candida glabrata, Candida krusei, SPC	24	400-035-C-MAX
	Syphilis	Treponema pallidum, SPC	24	400-038-C-MAX
Urinary Tract/Wound	Urinary Tract EKEGAS	Escherichia coli, Klebsiella spp., Enterococcus faecalis, Streptococcus pyogenes, SPC	24	400-028-C-MAX

CE-IVD Products (400-series part numbers) are not for sale in the US.
Controls and other info available at <https://www.biogx.com/syndromic-multiplexes/>

FDA-EUA

Syndromic Group	Product Name	Targets	No. of Reactions per kit	Part No.
COVID-19	Xfree COVID-19 Direct RT-PCR - Extraction-Free - Multi-Platform	SARS-CoV-2 N1, RNase P, IAC Extraction-Free Direct PCR - Multi-Platform	104 (direct)/160 (extracted) 52 (direct pooled)	500-003-XMP

Controls and other info available at <https://www.biogx.com/syndromic-multiplexes/>

BioGX US

1500 1st Avenue North, Birmingham,
Alabama, 35203, United States

Email: info@biogx.com

Phone: +1.205.250.8055

Website: www.biogx.com

BioGX BV

Science Park 408. 1098 XH
Amsterdam, The Netherlands

Email: eu@biogx.com

Phone: +31.20.893.426

Website: www.biogx.com



“LMP”, “HMP”, “XMP”, “XMM” and “PNP” are compatible with Applied Biosystems QuantStudio™ 5, Applied Biosystems 7500 Fast Dx, BioRad CFX96 Touch™, Bio-Rad CFX384 Touch™, Bio Molecular Systems Mic Instrument, BioGX pixl.16 or similar depending on the optical channel configuration of each product.

MAX’ product codes are compatible with BD MAX™ System from Becton Dickinson.

BD MAX trademark is the property of Becton Dickinson.

Applied Biosystems QuantStudio 5 trademark is the property of ThermoFisher Scientific.

Bio-Rad CFX96 Touch and Bio-Rad CFX384 Touch trademarks are the property of Bio-Rad Laboratories.