




## RAPID REAL-TIME PCR DETECTION OF MEASLES VIRUS

The resurgence of measles cases has underlined the need for fast and reliable testing solutions. BioGX Measles real-time PCR assay offers laboratories simplified testing solutions for ThermoFisher-ABI, Bio-Rad, and BD MAX™ platforms.

### BIOGX MEASLES ASSAY



- ✓ **Gold Standard:** Real-time PCR offers a rapid, non-invasive alternative to serology or viral culture.
- ✓ **Scalable Testing Format:** Enables low to high throughput screening of samples.
- ✓ **Comprehensive Screen:** Assay design provides an inclusive screen of the measles vaccine genotype and all measles circulating genotypes.
- ✓ **Flexible Deployment:** Designed for simplified integration into laboratories equipped with various real-time PCR platforms.



### BioGX Sample-Ready™ Real-Time PCR assay for the detection of Measles virus RNA (REF 450-114-LMP)\*

- ✓ **SPEED:** Sample-to-Result analysis in less than 2 hours.
- ✓ **PRECISION:** Multiplex PCR detection of Measles, RNase P, and RNA Internal Amplification Control. Assay design targets measles N-gene.
- ✓ **VERSATILITY:** Provides compatibility with ThermoFisher-ABI, QuantStudio™ Series, Bio-Rad CFX Series, and BD MAX™ system.
- ✓ **EFFICIENCY:** Simplified WorkFlow.

Contact Us



Email  
sales@biogx.com



Website  
www.biogx.com