



Amebic Infections

The BioGX Research Use Only real-time PCR assays for use with cerebral spinal fluid and homogenized corneal samples support rapid detection of protozoan pathogens.

Compatible sample types include:

- Cerebral Spinal Fluid and homogenized corneal samples

Product Specifications:

- REF 450-095-C-MAX
- Assay available in lyophilized format
- Compatible with BD MAX™
- Incorporation of sample processing control
- Positive Control available from BioGX (sold separately)
- Long-Term shelf-life stability

CONTACT US

✉ sales@biogx.com

🌐 www.biogx.com

Rev 1.0 Eff. Date 05JUNE2026

