


Candida auris

Sample-Ready™ Reagents for Real-time PCR 

The BioGX Research Use Only real-time PCR assay supports rapid detection of *Candida auris*. The assay is available in formats that support a variety of PCR instruments including BD MAX™, Bio-Rad CFX™, QuantStudio™, and pxiI™ real-time PCR platform.



Compatible sample types include:

- Urine (neat or boric acid preserved)
- Swab samples collected in Copan Universal Transport Media (UTM®) or Copan ESwab™ 
- Pretreated sputum samples
- Pretreated bronchoalveolar lavage (BAL) samples

Product Specifications:

- REF 450-043-Series
- BD MAX™: 2 hours for 24 samples
- PCR run after nucleic acid extraction on ABI or Bio-Rad platform: approx. 2hrs for 96 samples
- Incorporation of sample processing control
- Positive Control available from BioGX (sold separately)
- Long-Term shelf-life stability



BioGX Menu of *Candida auris* Assays

- *Candida auris* (450-043-Series)
- *Candida glabrata*, *Candida auris*, *Candida tropicalis* (450-096-LMP)



CONTACT US

 sales@biogx.com

 www.biogx.com

Rev 1.0 Eff. Date 05JUNE2026

