



Measles Virus

The BioGX Research Use Only real-time PCR assay supports rapid detection of RNA from Measles virus in upper respiratory tract samples. The assays are available in formats that support a variety of PCR instruments including BD MAX™, Bio-Rad CFX™, QuantStudio™, and pixl™ real-time PCR platform.



Compatible sample types include:

- Oropharyngeal Swab collections Copan Universal Transport Media (UTM®), Copan ESwab™, Viral Transport Media (VTM) , saline (0.85% NaCl), Phosphate Buffer Saline (PBS)



Product Specifications:

- REF 450-114-LMP
- Assay available in lyophilized format
- Multi-platform compatibility with Quantstudio 5™, ABI 7500 Fast Dx, Bio-Rad CFX96 Touch™, Bio-Rad CFX384 Touch™, pixl™ Real-Time PCR Platform, BD MAX™
- Incorporation of sample processing control
- Positive Control available from BioGX (sold separately)
- Long-Term shelf-life stability

CONTACT US

✉ sales@biogx.com

🌐 www.biogx.com

Rev 1.0 Eff. Date 05JUNE2026

